

# N736 Nitrogen

## Specification Sheet



<b>Instrument Range*</b>	0.0005 to 30 mg (0.5 ppm to 3.0% for a 1 g sample)	
<b>Precision**</b>	0.00025 mg (0.25 ppm) or 0.5% RSD, whichever is greater	
<b>Calibration</b>	Standards (single or multi-point); manual; gas dose <sup>†</sup>	
<b>Analysis Time</b>	He: 100 seconds; Ar: 130 seconds (including outgas, purge, and analysis delay)	
Cycle Time:	He: 180 seconds; Ar: 210 seconds (nominal)	
<b>Sample Size</b>	1 g (nominal)	
<b>Detection Method</b>	Thermal Conductivity	
<b>Chemical Reagents</b>	<ul style="list-style-type: none"><li>Anhydrous Magnesium Perchlorate (MgClO<sub>4</sub>)</li><li>Sodium Hydroxide on an Inert Base</li></ul>	<ul style="list-style-type: none"><li>Rare Earth Copper Oxide</li><li>Oxygen/Moisture Indicating Tube<sup>‡</sup></li><li>Copper<sup>‡</sup></li></ul>
<b>Gas Requirements</b>		
Carrier:	He: (99.99% pure), 22 psi (1.5 bar) ±5%	Ar: (99.999%), 22 psi (1.5 bar) ±5%
Pneumatic:	Compressed Air, 40 psi (2.8 bar) ±10%, source must be oil and water free	
<b>Gases Optional</b>		
Gas Dose <sup>†</sup> :	Nitrogen, 99.99% pure, 20 psi (1.4 bar) ±10%	
<b>Gas Flow Rates</b>		
Carrier:	480 cc/minute	
Pneumatic:	280 cc/analysis	
<b>Furnace</b>	Impulse furnace with current and power control 7500 Watts maximum, liquid cooled	
<b>Coolant</b>	3.2 L LECO Coolant	
<b>Operating Conditions</b>		
Temperature	15 to 35°C (59 to 95°F)	
Rel. Humidity	20 to 80%, non-condensing	
<b>Dimensions<sup>††</sup></b>		
Height:	36 in. (91.5 cm) nominal; 39.25 in. (100 cm) with load head cover lift engaged	
Width:	27.75 in. (71 cm)	
Depth:	30 in. (76 cm) without monitor; 31.5 in. (80 cm) with attached touch-screen monitor	
<b>Electrical Power</b>	230 V~ (+10/-15%; at max load); 50A, 50/60 Hz, Single Phase; 12,500 BTU/hr <sup>‡</sup>	
<b>Weight (approximate)</b>	Analyzer: 380 lb. (172 kg) without touch-screen monitor	

### Part Numbers

N736-XXXX Nitrogen Determinator with Windows®-based software and external PC

### Options

NOTE: Multiple configurations of options are available. Please contact your local LECO Sales Engineer for more details.

- Optional mounted touch-screen monitor package (M)
- Optional autocleaner package (H)
- Optional performance package (P)
- Optional dual cooling upgrade package (D)

\*Use the following formula to calculate element concentration:

$$\% \text{ element concentration} = ((\text{absolute element mass in mg})/(\text{sample mass in mg})) * 100$$

\*\*One sigma, conformance tested by gas dose analysis.

†Average output based on nominal operating parameters.

††Allow for a 6 in. (15 cm) minimum access area around all sides.

‡Optional.

**LECO**<sup>®</sup>  
Delivering the Right Results

## Theory of Operation


The N736 Nitrogen system is designed for measurement of nitrogen content of steel and other inorganic materials. The instrument features custom MS Windows®-based software designed specifically for touch operation.

A pre-weighed sample is placed in a graphite crucible which is heated in an impulse furnace to release analyte gases. Oxygen present in the sample reacts with the graphite crucible to form CO and CO<sub>2</sub>. An inert gas carrier, typically helium, sweeps the liberated gases out of the furnace and through a Mass Flow Controller. The gas then flows through a heated reagent, where the CO is oxidized to form CO<sub>2</sub>, and H<sub>2</sub> is oxidized to form H<sub>2</sub>O. CO<sub>2</sub> and H<sub>2</sub>O are then scrubbed out of the carrier gas stream. A Thermal Conductivity (TC) detector is used to detect the remaining nitrogen.

TC detection takes advantage of the difference in thermal conductivity between carrier and analyte gases. Resistive TC filaments are placed in a flowing stream of carrier gas and heated by a bridge circuit. As analyte gas is introduced into the carrier stream, the rate at which heat transfers from the filaments will change producing a measurable deflection in the bridge circuit.

The concentration of an unknown sample is determined relative to calibration standards. To reduce interferences from instrument drift, reference measurements of pure carrier gas are made prior to each analysis.

## Flow Diagram



Specifications and part numbers may change.  
Consult LECO for latest information.  
Windows® is a registered trademark of Microsoft Corporation.

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977  
info@leco.com • www.leco.com • ISO-9001:2008 HQ-Q-994 • LECO is a registered trademark of LECO Corporation.

LECO Corporation

**LECO**  
Delivering the Right Results