

O736 Oxygen

Specification Sheet



Instrument Range*	0.0005 to 2 mg (0.5 ppm to 0.2% for a 1 g sample)	
Precision**	0.00025 mg (0.25 ppm) or 0.5% RSD, whichever is greater	
Calibration	Standards (single or multi-point); manual	
Analysis Time	He: 85 seconds; Ar: 95 seconds (including outgas, purge, and analysis delay)	
Cycle Time:	He: 180 seconds; Ar: 190 seconds (nominal)	
Sample Size	1 g (nominal)	
Detection Method	Non-Dispersive Infrared Absorption	
Chemical Reagents	<ul style="list-style-type: none">• Anhydrous Magnesium Perchlorate (MgClO₄)• Sodium Hydroxide on an Inert Base• Rare Earth Copper Oxide, Copper Turnings• Oxygen/Moisture Indicating Tube[†]• Copper[†]	
Gas Requirements		
Carrier:	He: (99.99%), 22 psi (1.5 bar) ±5%	Ar: (99.999%), 22 psi (1.5 bar) ±5%
Pneumatic:	Compressed Air, 40 psi (2.8 bar) ±10%, source must be oil and water free	
Gases Optional		
Gas Dose [‡] :	Carbon Dioxide (CO ₂), 99.99% pure, 20 psi (1.4 bar/±10%)	
Gas Flow Rates		
Carrier:	450 cc/minute	
Pneumatic:	280 cc/analysis	
Furnace	Impulse furnace with current and power control 7500 Watts maximum, liquid cooled	
Coolant	3.2 L LECO Coolant	
Operating Conditions		
Temperature	15 to 35°C (59 to 95°F)	
Rel. Humidity	20 to 80%, non-condensing	
Dimensions^{††}		
Height:	36 in. (91.5 cm) nominal; 39.25 in. (100 cm) with load head cover lift engaged	
Width:	27.75 in. (71 cm)	
Depth:	30 in. (76 cm) without monitor; 31.5 in. (80 cm) with attached touch-screen monitor	
Electrical Power	230 V~ (+10/-15%; at max load); 50A, 50/60 Hz, Single Phase; 12,500 BTU/hr [†]	
Weight (approximate)	Analyzer: 380 lb. (172 kg) without touch-screen monitor	

Part Numbers

O736-XXXXC Oxygen Determinator with Windows®-based software and external PC

Options

NOTE: Multiple configurations of options are available. Please contact your local LECO Sales Engineer for more details.

- Optional mounted touch-screen monitor package (M)
- Optional autocleaner package (H)
- Optional performance package (P)
- Optional dual cooling upgrade package (D)

*Use the following formula to calculate element concentration:

$$\% \text{ element concentration} = ((\text{absolute element mass in mg})/(\text{sample mass in mg})) * 100$$

**One sigma, conformance tested by gas dose analysis.

[†]Average output based on nominal operating parameters.

^{††}Allow for a 6 in. (15 cm) minimum access area around all sides.

^{†††}Optional (with performance package P).

LECO[®]
Delivering the Right Results

Theory of Operation


The O736 Oxygen system is designed for measurement of oxygen content of steel and other inorganic materials. The instrument features custom MS Windows-based software designed specifically for touch operation.

A pre-weighed sample is placed in a graphite crucible which is heated in an impulse furnace to release analyte gases. Oxygen present in the sample reacts with the graphite crucible to form CO and CO₂. An inert gas carrier, typically argon or helium, sweeps the liberated gases out of the furnace and through a Mass Flow Controller. The gas flows through a heated reagent, where the CO is oxidized to form CO₂, and H₂ is oxidized to form H₂O. Oxygen is detected as CO₂ using a non-dispersive infrared (NDIR) cell. CO₂ and H₂O are then scrubbed out of the carrier gas stream.

NDIR cells are based on the principle that analyte gas molecules absorb infrared (IR) energy at unique wavelengths within the IR spectrum. Incident IR energy at these wavelengths is absorbed as the gases pass through the IR absorption cells.

The concentration of an unknown sample is determined relative to calibration standards. To reduce interferences from instrument drift, reference measurements of pure carrier gas are made prior to each analysis.

Flow Diagram



Specifications and part numbers may change.
Consult LECO for latest information.
Windows® is a registered trademark of Microsoft Corporation.

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977
info@leco.com • www.leco.com • ISO-9001:2008 HQ-Q-994 • LECO is a registered trademark of LECO Corporation.

LECO Corporation

LECO[®]
Delivering the Right Results